




Innovating Technologies for Life

CrocClamp™

Non-Reopening Tubing Clamp

FEATURES

- ▶ *Ergonomic design*
- ▶ *Easy to close*
- ▶ *Locking closure mechanism*
- ▶ *No accidental clamp release*
- ▶ *Permanent seal*



CrocClamp™



Non-Reopening Tubing Clamp

Instructions for use




Place fingers along the top and bottom of the CrocClamp so that they extend beyond the edge of the clamp.



Squeeze together the top and bottom of the CrocClamp to ensure full closure of the clamp. This will provide occlusion of the tubing.



CrocClamp in locked position.



For information please contact



Innovating Technologies for Life

ITL Limited (Head Office)

t: +(61) 2 6228 3300
f: +(61) 2 6228 3377
sales@itl-limited.com

www.itl-limited.com

ITL Australia

t: +(61) 3 8773 3000
f: +(61) 3 8773 3091
sales@itlhealthcare.com

ITL North America and Europe

t: +(1) 703 435 6700
f: +(1) 703 435 6717
sales@itlus.com

ITL Asia Pacific

t: +(60) 3 8023 0486
f: +(60) 3 8023 0487
sales@itlcorporation.com